

## PERFORMANCE ♦ EFFICIENCY ♦ VALUE

MEDTOX<sup>®</sup> screening devices provide accurate, easy-to-use, and cost-effective testing for drugs-of-abuse. Since 1984 MEDTOX Diagnostics, Inc. has been developing innovative ways to assist hospitals and occupational health facilities with their toxicology and drugs-of-abuse testing. MEDTOX is pleased to bring you the PROFILE II ER with visual read. Simply add specimen and read at 5 - 15 minutes.



### MEDTOX Screening Levels:

AMP	1000 ng/ml
BAR	200 ng/ml
BZO	300 ng/ml
COC	300ng/ml
<b>MAMP</b>	<b>1000 ng/ml</b>
<b>MDMA</b>	<b>1500 ng/ml</b>
MTD	300 ng/ml
OPI	300 ng/ml
<b>OXY</b>	<b>100 ng/ml</b>
PCP	25 ng/ml
PPX	300 ng/ml
TCA	300 ng/ml
THC	50 ng/ml

### Ordering Information for PROFILE-II Devices

Cat. No.	Description	U/M
601622	PROFILE-II ER 7-Panel Device (AMP, BAR, BZO, COC, OPI, PCP, THC)	25/BX
601591	PROFILE-II ER 9-Panel Device (AMP, BAR, BZO, COC, MTD, OPI, PCP, TCA, THC)	25/BX
601674	PROFILE-II ER 12-Panel Device* (AMP, BAR, BZO, COC, MAMP/MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC) <i>*Note that this device includes Oxycodone</i>	25/BX
101830	MEDTOX Positive Control (OPI 300mg/ml), 5 ml bottle	Each
101183	MEDTOX Negative Control 5ml bottle	Each