

PERFORMANCE ♦ EFFICIENCY ♦ VALUE

PROFILE[®] V MEDTOXScan[®] Drugs of Abuse Test System

The next generation in hospital laboratory drug screening

- Easy installation and operation
- One-step testing process
- Results in ten minutes
- Lower clinical cut-off levels
- Stores up to 1000 results.
- Eliminates operator result interpretation
- Can test up to 12 drugs, including Methamphetamine and Oxycodone
- Ability to report into hospital LIS



PROFILE[®]-V MEDTOXScan[®] Drugs of Abuse Test System

MEDTOX[®] screening devices and reader provide accurate, easy-to-use, and cost-effective testing for drugs-of-abuse. Since 1984 MEDTOX Diagnostics, Inc. has been developing innovative ways to assist hospitals and occupational health facilities with their toxicology and drugs-of-abuse testing. MEDTOX is pleased to bring you the PROFILE[®]-V MEDTOXScan[®] Drugs of Abuse Test System.

MEDTOXScan Reader Features:

- Barcode scanner for accurate entry of patient identification, specimen ID and lot number
- Provides a digital display of test results
- Printer provides paper copy of results
- Provides operator identifications
- Automated device detection
- Ability to report into hospital LIS via middleware

MEDTOXScan uses the MEDTOX PROFILE[®]-V devices. The devices are currently available in a 12, 9 and 7 panel.

PROFILE[®]-V Screening Levels:

| | |
|------|-----------|
| AMP | 500 ng/ml |
| BAR | 200 ng/ml |
| BZO | 150 ng/ml |
| COC | 150 ng/ml |
| MAMP | 500 ng/ml |
| MTD | 200 ng/ml |
| OPI | 100 ng/ml |
| OXY | 100 ng/ml |
| PCP | 25 ng/ml |
| PPX | 300 ng/ml |
| TCA | 300 ng/ml |
| THC | 50 ng/ml |

Ordering Information for MEDTOXScan

| Cat. No. | Description | U/M |
|----------|---|----------|
| 833073 | MEDTOXScan Device Reader | Each |
| 833063 | MEDTOXScan Printer | Each |
| 101885 | Printer Paper (Replacement) | Each |
| 833061 | MEDTOXScan Barcode Reader | Each |
| 101885 | Paper Rolls for MEDTOXScan | Each |
| 833075 | MEDTOXScan QC Test Devices, Positive and Negative | Each set |
| 833076 | Cleaning Cassette | Each |

Ordering Information for PROFILE-V Devices

| Cat. No. | Description | U/M |
|----------|---|-------|
| SP604020 | PROFILE-V 7-Panel Device (AMP, BAR, BZO, COC, OPI, PCP, THC) | 25/BX |
| SP604019 | PROFILE-V 9-Panel Device (AMP, BAR, BZO, COC, MTD, OPI, PCP, TCA, THC) | 25/BX |
| SP604018 | PROFILE-V 12-Panel Device (AMP, BAR, BZO, COC, MAMP, MTD, OPI, OXY, PCP, PPX, TCA, THC) | 25/BX |
| SP102081 | PROFILE-V MEDTOX Scan Positive Control (OPI 300mg/ml), 5 ml bottle | Each |
| SP101183 | MEDTOX Negative Control 5ml bottle | Each |

MEDTOX controls must be used with MEDTOX devices.

