

**PERFORMANCE ♦ EFFICIENCY ♦ VALUE**

## **PROFILE<sup>®</sup> V MEDTOXScan<sup>®</sup> Drugs of Abuse Test System**

**The next generation in hospital laboratory drug screening**

- Easy installation and operation
- One-step testing process
- Results in ten minutes
- Lower clinical cut-off levels
- Stores up to 1000 results.
- Eliminates operator result interpretation
- Can test up to 12 drugs, including Methamphetamine and Oxycodone
- Ability to report into hospital LIS



**PROFILE<sup>®</sup>-V MEDTOXScan<sup>®</sup> Drugs of Abuse Test System**

MEDTOX<sup>®</sup> screening devices and reader provide accurate, easy-to-use, and cost-effective testing for drugs-of-abuse. Since 1984 MEDTOX Diagnostics, Inc. has been developing innovative ways to assist hospitals and occupational health facilities with their toxicology and drugs-of-abuse testing. MEDTOX is pleased to bring you the PROFILE<sup>®</sup>-V MEDTOXScan<sup>®</sup> Drugs of Abuse Test System.

**MEDTOXScan Reader Features:**

- Barcode scanner for accurate entry of patient identification, specimen ID and lot number
- Provides a digital display of test results
- Printer provides paper copy of results
- Provides operator identifications
- Automated device detection
- Ability to report into hospital LIS via middleware

MEDTOXScan uses the MEDTOX PROFILE<sup>®</sup>-V devices. The devices are currently available in a 12, 9 and 7 panel.

**PROFILE<sup>®</sup>-V Screening Levels:**

AMP	500 ng/ml
BAR	200 ng/ml
BZO	150 ng/ml
COC	150 ng/ml
MAMP	500 ng/ml
MTD	200 ng/ml
OPI	100 ng/ml
OXY	100 ng/ml
PCP	25 ng/ml
PPX	300 ng/ml
TCA	300 ng/ml
THC	50 ng/ml

**Ordering Information for MEDTOXScan**

Cat. No.	Description	U/M
833073	MEDTOXScan Device Reader	Each
833063	MEDTOXScan Printer	Each
101885	Printer Paper (Replacement)	Each
833061	MEDTOXScan Barcode Reader	Each
101885	Paper Rolls for MEDTOXScan	Each
833075	MEDTOXScan QC Test Devices, Positive and Negative	Each set
833076	Cleaning Cassette	Each

**Ordering Information for PROFILE-V Devices**

Cat. No.	Description	U/M
SP604020	PROFILE-V 7-Panel Device (AMP, BAR, BZO, COC, OPI, PCP, THC)	25/BX
SP604019	PROFILE-V 9-Panel Device (AMP, BAR, BZO, COC, MTD, OPI, PCP, TCA, THC)	25/BX
SP604018	PROFILE-V 12-Panel Device (AMP, BAR, BZO, COC, MAMP, MTD, OPI, OXY, PCP, PPX, TCA, THC)	25/BX
SP102081	PROFILE-V MEDTOX Scan Positive Control (OPI 300mg/ml), 5 ml bottle	Each
SP101183	MEDTOX Negative Control 5ml bottle	Each

*MEDTOX controls must be used with MEDTOX devices.*

